

SPONGE, NOVAPOLY

Product Specification

Product Name:	Sponge, NovaPoly
Manufacturer Part Number(s):	SP6-33, SP6-363, SP6-46, SP6-745
Description:	A smooth knit polyester cover over a urethane foam sponge for cleaning small parts, equipment, and other lab applications. Available in a variety of shapes and sizes. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Outer Layer	100% smooth knit polyester
Interior	Urethane foam
Dimensions:	
Width x Length x Height	SP6-33: 3 in x 3 in x 0.5 in (7.6 cm x 7.6 cm x 1.2 cm) SP6-363: 3 in x 6 in x 3 in (7.6 cm x 15 cm x 7.6 cm) SP6-46: 4 in x 6 in x 0.5 in (10 cm x 15 cm x 1.2 cm) SP6-745: 7 in x 4.5 in x 0.5 in (17.8 cm x 11.4 cm x 1.2 cm)
For Use with: (as needed)	N/A
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions
Packaging:	Individually bagged; 24 sponges per case
Traceability:	Lot numbers printed on each bag and each case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website