

## Silsat Wipes, Irradiated

### Product Specification

Product Name:	Silsat Wipes, Irradiated
Manufacturer Part Numbers:	SSW6-01IR
Description:	An irradiated knit polyester wipe, pre-saturated with 50% medical grade silicone and 50% IPA for use in aseptic processing and support areas to lubricate gaskets and seals in isolators, refrigerators, and cabinets. Packaged 8 wipes per resealable pack. Manufactured at Micronova's Torrance, CA Facility. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Wipe fabric: 100% knit polyester Saturation: 50% silicone, 50% alcohol
Dimensions:	6 inches x 4 inches x 0.5 inches (15.2 cm x 10 cm x 1.2 cm)
For Use with: (as needed)	Aseptic processing and support areas equipment
Autoclavability:	Not autoclavable
Chemical Compatibility:	Good compatibility with common disinfectants and cleaning solutions
Shelf life:	18 months (inclusive of manufacturing date)
Packaging:	8 wipes pre resealable pack, 24 packs per poly lined case. Each pack is double bagged
Traceability:	Lot numbers and part numbers printed on each pack and outer case. Irradiation indicator placed on each outer bag and each case.
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certifications:	Certificate of Conformance and Certificate of Irradiation Processing shipped with each order and available on website. Sterility test results available upon request.