

CLEANROOM POLYETHYLENE TAPE

Product Specification

Product Name:	Cleanroom Polyethylene Tape
Manufacturer Part Number(s):	PCX-125DB, PCX-15TR, PCX-1BL, PCX-1DB, PCX-1DG, PCX-1OR, PCX-1RE, PCX-1TR, PCX-1WH, PCX-1YE, PCX-2BL, PCX-2DB, PCX-2DG, PCX-2OR, PCX-2RE, PCX-2TR, PCX-2WH, PCX-2YE, PCX-3OR, PCX-34BL, PCX-34DB, PCX-34DG, PCX-34OR, PCX-34RE, PCX-34TR, PCX-34WH, PCX-34YE, PCX-4WH
Description:	<p>A low-density polyethylene tape that is low in halogens and leachable chlorides. Available in a variety of colors and widths, the tape can be used for sealing, color-coding, labeling, batching and food processing applications.</p> <p>Manufactured at Micronova's Torrance, CA facility.</p> <p>Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.</p>
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Film	Polyethylene
Adhesive	Acrylic
Core	Plastic
Color	BL: Black, DB: Dark Blue, DG: Dark Green, OR: Orange RE: Red, TR: Transparent, WH: White, YE: Yellow
Dimensions:	3 in (7.6 cm) 36 yards (33 m) PCX-34: ¾ in (1.9 cm) PCX-125: 1.25 in (3.1 cm) PCX-1: 1 in (2.5 cm) PCX-15: 1.5 in (3.8 cm) PCX-2: 2 in (5 cm) PCX-3: 3 in (7.6 cm) PCX-4: 4 in (10 cm)
For Use with: (as needed)	P-52, P-52W Tape Dispensers
Autoclavability:	Not autoclavable
Chemical Compatibility:	Good chemical compatibility with common solutions
Packaging:	Individually bagged PCX-34: 16 rolls per case PCX-1: 36 rolls per case PCX-15: 24 rolls per case PCX-2 and PCX-125: 24 rolls per case PCX-3 and PCX-4: 12 rolls per case
Traceability:	Part number and lot number printed on each bag and case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website

