



SLIPCOVER MICROFIBER

Product Specification

Product Name:	SlipCover Microfiber
Manufacturer Part Numbers:	SMCAA-14, SMCAA-18
Description:	SlipCover for use with SlimLine, SnapMop or PadMop Series. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001:2008 Facility
Physical Data:	100% Microfiber
Measurement:	
SMCAA-14:	5" x 14"
SMCAA-18:	5" x 18"
Packaging:	SlipCovers are individually bagged and packaged 50 SlipCovers per poly lined case.
Traceability:	Lot numbers printed on each outer bag and each case.