



ISOLATOR TOOL COVER - POLYSORB

Product Specification

Product Name:	Isolator Tool Cover PolySorb
Manufacturer Part Numbers:	ITC11-9NS, ITC11-14NS, ITC11-18NS
Description:	<p>A two ply textured polyester slip cover for use with the QDIT and QPIT series mop adapters. Available in three sizes 9", 14" and 18" Manufactured at Micronova's Torrance, CA facility.</p> <p>Designed for isolators, bio-hazard and containment hoods</p>
Manufacturing Standards:	ISO 9001:2008 Facility
Physical Data:	100% Polyester with a textured 'waffle weave' for increased absorbency
Packaging:	Each mop head single bagged and packaged 50 mops per poly lined case.
Traceability:	Lot numbers printed on each bag and each case