



ISOLATOR TOOL COVER – POLYSORB IRRADIATED

Product Specification

| | |
|----------------------------|--|
| Product Name: | Isolator Tool Cover PolySorb Irradiated |
| Manufacturer Part Numbers: | ITC11-9NSIR, ITC11-14NSIR, ITC11-18NSIR |
| Description: | <p>A two ply textured polyester slip cover for use with the QDIT and QPIT series mop adapters. Available in three sizes 9", 14" and 18" Manufactured at Micronova's Torrance, CA facility.</p> <p>Designed for isolators, bio-hazard and containment hoods</p> |
| Manufacturing Standards: | ISO 9001:2008 Facility |
| Physical Data: | 100% Polyester with a textured 'waffle weave' for increased absorbency |
| Packaging: | <p>Each mop head packaged in Class 100 cleanroom. Double bagged. Irradiation indicator placed on outer bag. 50 mops are then packed into a poly lined case. Irradiation indicator placed on case.</p> |
| Traceability: | Lot numbers printed on each outer bag and each case |
| Irradiation Specification: | Minimum irradiation dose 25 kGy |
| Certificates: | Certificate of irradiation shipped with each order. |