



## ISOLATOR TOOL

### Product Specification

Product Name:	Isolator Tool
Manufacturer Part Numbers:	QDIT-09, QDIT-14, QDIT-18
Description:	An electropolished stainless steel frame with a High-density polyethylene (HDPE) swivel joint. For use with ITC series of polyester slip cover mops. Available in three sizes 9", 14" and 18" Manufactured at Micronova's Torrance, CA Facility Designed for isolators, bio-hazard and containment hoods
Manufacturing Standards:	ISO 9001:2008 Facility
Physical Data:	304 Stainless steel - electropolished Polypropylene connector unit Stainless steel and delrin fittings
Packaging:	Each Adapter Tool single bagged and packaged 6 units per case