



SWABHOLDER ADAPTER

Product Specification

Product Name:	SwabHolder Adapter
Manufacturer Part Numbers:	QDSWT-01
Description:	Electropolished stainless steel adapter securely holds environmental monitoring swabs. Quick Disconnect attaches to SSU series handles for extended reach. Manufactured at Micronova's Torrance, CA facility Designed for cleanroom use
Manufacturing Standards:	ISO 9001:2008 facility
Materials:	304 Stainless Steel – Electropolished with Polypropylene and Delrin Fittings
Packaging:	Each adapter single bagged and packaged 6 per case.