



PHARMAPLY MOP IRRADIATED

Product Specification

Product Name:	PharmaPly Mop Irradiated
Manufacturer Part Numbers:	PPM4-MM15IR
Description:	A poly/rayon string mop for use with the QDMH-1 or MH-1 series mop adapters. Designed for Controlled Environments and Cleanroom Support applications
Physical Data:	Polyester/Rayon
Packaging:	Each mop head is packaged in a cleanroom. Double bagged in groups of 4.
Irradiation case.	indicator placed on outer bag and outer case. 24 mops packed in a poly lined case.
Traceability:	Lot numbers printed on each outer bag and each case
Irradiation Specification:	Minimum irradiation dose 40 kGy
Certificates:	Certificate of irradiation shipped with each order