



ISOLATOR TOOL COVER MEGATEX IRRADIATED

Product Specification

Product Name:	Isolator Mop Cover MegaTex Irradiated
Manufacturer Part Numbers:	ITC620-9IR, ITC620-14IR, ITC620-18IR
Description:	<p>An elasticized two ply polyester slip cover with a strip of MegaTex for use with the QDIT/QPIT series mop adapters. MegaTex is a non-woven polyamide fabric with a textured thermo-bonded surface. Designed for isolators, bio-hazard and containment hoods</p> <p>Available in three sizes 9", 14" and 18"</p> <p>Manufactured at Micronova's Torrance, CA facility.</p>
Manufacturing Standards:	<p>ISO 9001:2008 Facility</p> <p>10⁻⁶ Sterility Assurance Level</p> <p>In-compliance with AAMI Guidelines for Sterilization</p>
Physical Data:	Polyester/Polyamide
Packaging:	<p>Each mop head packaged in a cleanroom. Double bagged. Irradiation indicator placed on outer bag. 50 mops packed into a poly lined case.</p> <p>Irradiation indicator placed on case.</p>
Traceability:	Lot numbers printed on each outer bag and each case
Irradiation Specification:	Minimum irradiation dose 40 kGy
Certificates:	<p>Certificate of irradiation shipped with each order.</p> <p>Validation paperwork available on request.</p>