

SilSat™ Sterile Lubricating Wipes



The Micronova SilSat™ Sterile Lubricating Wipes are designed for use in aseptic processing and support areas to lubricate gaskets and seals in isolators, refrigerators, and cabinets.

The irradiated knit polyester wipe is pre-saturated with 50% medical grade silicone and 50% IPA for a balanced lubrication, without over saturation.

- Sterile, pre-saturated wipes
- Designed to lubricate gaskets, seals and equipment
- Gamma irradiated for use in aseptic areas
- Processed for cleanrooms and controlled environments, ISO 5 to ISO 9
- Good compatibility with common disinfectants and cleaning solutions



SPECIFICATIONS

Materials of Construction:

50% IPA

Dimensions:

Wipe: 6 x 4 inch, 15.2 x 10 cm
Package Size : 7-1/2"(L) x 4 - 1/4"(W)

Packaging Configurations:

- 8 wipes per resealable pack, 24 packs per poly-lined case; each pack is double bagged
- 4 wipes pre resealable pack, 48 packs per poly-lined case; each pack is double bagged

Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO 9001 Certified Manufacturing Facility
- Minimum irradiation dose: 25 kGy
- Store at ambient temperature
- Avoid exposure to direct sunlight
- Shelf Life: 18 months [inclusive of manufacturing date]



Micronova Manufacturing
3431 W. Lomita Boulevard
Torrance, CA 90505-5010 USA

Tel: 310/784-6990
Fax: 310/784-6980

SalesInfo@micronova-mfg.com
www.micronova-mfg.com