

# LYOTOOL

## Product Specification

Product Name:	LyoTool
Manufacturer Part Number(s):	QDBP14-02, QDBP14-03, QDBP14-04, QDBP18-03
Description:	An adjustable dual plastic plate mop adapter holds SnapMops, PocketMops, or PadMops to effectively cleans top and bottom of shelves simultaneously. Attaches to stainless SSU series stainless steel handles. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Plates and Posts:	Acetal
Handle Connection:	Stainless Steel and Acetal
Dimensions:	
Plates:	QDBP14 series: 14 in x 4 in (350 mm x 100 mm) QDBP18 series: 18 in x 4 in (450 mm x 100 mm)
Height Range:	QDBP14-02: Adjusts 2 in to 3 in max height (114 mm to 152 mm) QDBP14-03, QDBP18-03: Adjusts 3 in to 4 ½ in max height (76 mm to 144 mm) QDBP14-04: Adjusts 4 in to 6 in max height (114 mm to 152 mm)
For Use with: (as needed)	SSU series handles and SnapMops, PocketMops or PadMops as covers
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Packaging:	Individually bagged
Traceability:	Part Number and Lot Number on each bag
Irradiation Standards:	Not irradiated
Certificates:	Certificate of Conformance available on request