

Yellow Biohazard Bag, Irradiated

Product Specification

Product Name:	Yellow Biohazard Bag, Irradiated
Manufacturer Part Numbers:	BL2-32YEIRBB.
Description:	Irradiated yellow biohazard bag LDPE, 2 mil, Star Bottom. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	LLDPE Yellow Bag with black printed BioHazard symbol
Dimensions:	
Thickness:	2 mil
Bag size:	32 x 36 in (80 x 90 cm)
For Use with:	BLSSH-1515
Autoclavability:	Not Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Packaging:	5 liners per double bag
Traceability:	Lot numbers printed on each package and each case.
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certifications:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.