

# LyoTool™

## Unique Dual Action Lyophilizer Cleaning



The new adjustable **LyoTool™** with dual **SnapMops™** effectively cleans top and bottom of lyophilizer or incubator shelves with one pass. Compression fit plates ensure even sanitization and consistent contact with plates as liquid leaves the mop surface.

Designed to clean lyophilizer shelves, the tool works well on a variety of equipment with shelves spaced from 3 to 6 inches apart. The LyoTool mop plates are adjustable to match the height of the lyophilizer shelves with a simple twist of the knob.



Shown with SnapMops™ installed ▶

- Two popular sizes for common shelf height, adjustable 3 to 4 1/2 inches or 4 1/2 to 6 inches.
- Available 14" and 18" width to match aseptic mops.
- Front Nose Guard to keep tool on the shelf.
- Spring loaded shelves for consistent pressure up and down.
- QD series for Quick Disconnect SSU series stainless steel handles.
- QP series for ACME threaded STA series plastic handles.
- Smooth NovaPoly™ SnapMop covers are ideal for dispensing IPA or other sanitizers.
- PolySorb™ or PolyMesh™ SnapMop covers are excellent for cleaning product residue from the shelves, as well as sanitizing.
- SnapMops are autoclavable and available gamma irradiated.

- QDBP14-03 LyoTool, 14 inch, 14 X 4 X 4-1/2"  
Adjusts 3 to 4-1/2" max. height (76 to 114mm)
- QDBP14-04 LyoTool, 14 inch, 14 X 4 X 6"  
Adjusts 4-1/2 to 6" max height (114 to 152mm)
- QDBP18-03 LyoTool, 18 inch, 18 X 4 X 4-1/2"  
Adjusts 3 to 4-1/2" max height (76 to 114mm)
- QDBP18-04 LyoTool, 18 inch, 18 X 4 X 6"  
Adjusts 4-1/2 to 6" max height (114 to 152mm)

### Materials of Construction:

Plates and Posts: Acetal  
Handle Connection: QD: Stainless and Acetal  
QP: Acetal

**Weight:** 2 pounds / .9 kilograms

Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO Certified Manufacturing Facility
- Autoclavable
- Good chemical resistance to disinfectants



Micronova Manufacturing, Inc.  
3431 W. Lomita Boulevard  
Torrance, CA 90505-5010 USA

Tel: 310/784-6990  
Fax: 310/784-6980

SalesInfo@micronova-mfg.com  
www.micronova-mfg.com