

AUTOCLAVE INDICATOR TAPE, LEAD FREE

Product Specification

Product Name: Autoclave Indicator Tape, Lead Free

Manufacturer Part Number(s): AUTOPC-1, AUTOPC-12, AUTOPC-15, AUTOPC-2, AUTOPC-25,

AUTOPC-3, AUTOPC-34

Description: A lead-free indicator tape for steam sterilization cycles. Green

indicator lines change to black at completion of the sterilization

cycle in both gravity discharge and vacuum autoclaves. Manufactured at Micronova's Torrance, CA facility.

Designed for cleanrooms and controlled environments, ISO 5 to

ISO 9.

Manufacturing Standards: ISO 9001 Certified Facility

RoHS compliant

Physical Data:

Materials of Construction:

Film Saturated crepe paper Adhesive Synthetic elastomer

Core Plastic Color Natural

Dimensions:

Core 3 in (7.6 cm) Length 60 yards (55 m)

Width AUTOPC-34: ¾ in (1.9 cm)

AUTOPC-12: ½ in (1.3 cm) AUTOPC-1: 1 in (2.5 cm) AUTOPC-15: 1.5 in (3.8 cm) AUTOPC-2: 2 in (5 cm) AUTOPC-25: 2.5 in (6.4 cm) AUTOPC-3: 3 in (7.6 cm)

For Use with: (as needed) P-52, P-52W Tape Dispensers

Autoclavability: Autoclavable

Chemical Compatibility: Good chemical compatibility with common solutions

Packaging: Individually bagged

AUTOPC-1: 36 rolls per case AUTOPC-12: 72 rolls per case

AUTOPC-15, AUTOPC-2: 24 rolls per case AUTOPC-25, AUTOPC-3: 12 rolls per case

AUTOPC-34: 48 rolls per case

Traceability: Part number and lot number printed on each bag and case

Irradiation Standards: Not Irradiated

Certificates: Certificate of Conformance available upon request and on

website

Phone: 310-784-6990 / micronova-mfg.com