

CLEANROOM VINYL TAPE, IRRADIATED

Product Specification

Product Name:	Cleanroom Vinyl Tape, Irradiated	
Manufacturer Part Numbers:	CR100-xxx IR, where "xxx" is the size and color, such as CR100-1WHIR for one inch white, or CR100-12WHIRM, for ½ inch white, small case. See sizes and case pack below.	
Description:	An ultra-processed, irradiated vinyl tape wound on a 3" plastic core and 36 yards in length. Available in a variety of widths. Designed for cleanrooms and controlled environments, ISO5 to ISO9.	
Manufacturing Standards:	ISO 9001 Certified Facility	
Physical Data:		
Materials of Construction:	Backing: White, pigmented and transparent vinyl film Adhesive: Rubber based pressure sensitive	
Dimensions:	Width: varies by part number Length: 36 yards	
For Use with: (as needed)	Good adhesion and removal from stainless steel, glass, and epoxy surfaces	
Autoclavability:	Autoclavable when applied to glass or stainless-steel surfaces	
Chemical Compatibility:	Good chemical compatibility with common solutions	
Shelf life:	24 months	
Packaging:	Individually double bagged	
	Part Number Suffix – Size, Color, IRM:	Part Number Suffix – Size, Color, IR
	3/8" = 32 rolls per case	3/8" = 96 rolls per case
	1/2" = 24 rolls per case	1/2" = 72 rolls per case
	3/4" = 16 rolls per case	3/4" = 48 rolls per case
		1" = 36 rolls per case
		2" = 24 rolls per case
		3" = 12 rolls per case
Traceability:	Lot numbers printed on each roll and ease case	
Irradiation Standards:	Minimum irradiation dose 25kGy	
Certifications:	Certificate of Irradiation and Certificate of Conformance shipped with each order	