

Bucket Binder Double Bucket Cart System with Stainless Steel Casters

Product Specification

Product Name:	Bucket Binder Double Bucket Cart System with Stainless Steel Casters
Manufacturer Part Number(s):	BB-DBL-DBYESS, BB-DBL-GRYESS
Description:	A modular double plastic bucket cart system with heavy duty stainless-steel caster that is easy to use, maneuver and store. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Bucket Binders:	Electropolished Stainless Steel and Polyketone
Buckets:	Polypropylene
Casters:	Stainless Steel
Wheels:	Polyolefin
Dimensions:	
Length:	22 in (56 cm)
Width:	22 in (56 cm)
Height:	14.5 in (37 cm)
For Use with: (as needed)	W-94, W-9, W-92, W-9B
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Shelf Life:	N/A
Packaging:	Individually Bagged
Traceability:	Lot number and part number printed on each bag and each case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request