

BUCKETLINER

Product Specification

Product Name:	BucketLiner
Manufacturer Part Number(s):	BKL7-3119, BKL-2424
Description:	LDPE/Nylon laminate liners from FDA compliant LDPE with a horizontal orientation to be used as an impermeable barrier for lining, buckets and containers. The liner can also be used to hold liquids going through the autoclave. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	Clear LDPE 75%, nylon 25% laminate
Dimensions: Thickness Open Width x Height Capacity	3 mil BKL7-3119: 31 in x 19 in (79 cm x 28 cm) BKL-2424: 24 in x 24 in (60 cm x 60 cm) 6 gallons, 22 liters
For Use with: (as needed)	Use BKL7-3119 with: B-7, B-7AC series buckets and BucketBinder Multi Bucket Systems Use BKL-2424 with C-2GRPW and C-2PW
Autoclavability:	Autoclavable - Will hold up to autoclaving with liquid contents
Chemical Compatibility:	Good compatibility with common cleaners and solutions
Shelf Life:	NA
Packaging:	Flat packed, 20 liners per bag, 200 liners per case
Traceability:	Lot numbers printed on each package and each case.
Irradiation Standards:	Not irradiated
Certificates:	Certificate of Conformance available on request