

POLYMITT™ CLEANING MITT WIPE, MICROFIBER

Product Specification

Product Name:	PolyMitt™ Cleaning Mitt Wipe, Microfiber
Manufacturer Part Number(s):	MW-AA, MW-AW, MWDT-AA11
Description:	A microfiber cleaning mitt designed to fit over a gloved hand for cleaning around carts, shelving, tanks, tubing and other hard to clean areas. Microfiber material is highly absorbent and effective wet or dry. Available in one standard size, with or without a thin urethane foam lining. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Outer Layer	80% polyester, 20% polyamide
Color	MW-AA, MWDT-AA11: Blue MW-AW: White
Inner Layer	MW-AA/AW: No foam MWDT-AA11: Thin urethane foam; both sides
Dimensions:	
Length x Width	12 in x 8 in (30.5 cm x 20.3 cm)
Thickness	MW-AA/AW: 0.25 in (0.6 cm) MWDT-AA11: 0.50 in (1.3 cm)
For Use with: (as needed)	Carts, shelving, tanks, tubing
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Shelf Life	N/A
Packaging:	Individually bagged; 24 mitts per poly lined case
Traceability:	Lot number and part number printed on each bag and each case
Irradiation Standards:	Not irradiated
Certificates:	Certificate of Conformance available upon request