



WETDRY WIPES

Product Specification

Product Name:	WetDry Wipes
Manufacturer Part Numbers:	MZSW40-910
Description:	Combo pack. 2 individually packaged Sontara Wipes – 1 wet/1 dry. Wet wipe saturated with custom blend of Lab & Glass detergent and IPA (80/20%). For lens and instrument cleaning.
Physical Data:	8" x 9" Sontara Wipe packaged in a 3" x 5" foil pack (x 2 perforated attachment).
Packaging:	200 packs (wet and dry combo) per poly lined case. Part number and identification (wet/dry) per combo pack. Part number on outer case
Traceability:	Lot numbers printed on each case