

# M-ZONE LOOPMOP, POLYTUBE, IRRADIATED

## Product Specification

Product Name:	M-Zone LoopMop, PolyTube, Irradiated
Manufacturer Part Number(s):	MZLM5-12IR
Description:	An irradiated, 100% tubular knit polyester looped string mop with no exposed raw edges to ensure low particulation for applying disinfectant in pharmaceutical, aseptic processing and medical device applications. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments, ISO 5 to ISO 9.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Strands	100% polyester tubular knit
Band	Polyester
Dimensions:	
Length	15 in (38 cm)
Weight	12 oz
For Use with: (as needed)	MH-1 or MHC-1 fiberglass handles, or QDMH-1 MicroMop frame with SSU series handle
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Shelf Life:	60 months
Packaging:	Individually double bagged, 24 mop heads in a lined case
Traceability:	Lot numbers printed on bag and each case
Irradiation Standards:	Minimum irradiation dose 25 kGy
Certificates:	Certificate of Irradiation Processing and Certificate of Conformance shipped with each order and available on website.