

PADMOP, NOVAPOLY

Product Specification

Product Name:	PadMop, NovaPoly
Manufacturer Part Number(s):	SLP6-9, SLP6-14, SLP6-18
Description:	<p>A brick style, smooth polyester covered sponge pad mop, with a foam layer both above and below the mop frame. Designed to apply the optimum amount of disinfectant in pharmaceutical production and aseptic areas. Manufactured at Micronova's Torrance, CA facility. Designed for cleanrooms and controlled environments, ISO 5 to ISO 9.</p>
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Outer layer	100% smooth knit polyester
Interior	Urethane foam
Dimensions:	
Total Thickness	1.5 in (3.8 cm)
Length x Width	SLP6-9: 9 in x 5 in (22.8 cm x 12.7 cm) SLP6-14: 14 in x 5 in (35.5 cm x 12.7 cm) SLP6-18: 18 in x 5 in (45.7 cm x 12.7 cm)
For Use with: (as needed)	QPSL or QSSL frame with STA series handle or QDSL with SSU series handle
Autoclavability:	Autoclavable
Chemical Compatibility:	Good compatibility with standard disinfectants and cleaning solutions
Shelf Life:	NA
Packaging:	Individually bagged; 12 mop heads in a lined case
Traceability:	Lot numbers printed on each bag and each case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request