

TMOPT™, NOVAPOLY

Product Specification

Product Name:	TMop™, NovaPoly
Manufacturer Part Number(s):	T6-10, T6-14, T6-18
Description:	A smooth knit polyester mop cover with a foam interior and stainless steel fasteners to be used with the QDT, QPT and QPWT series mop adapters. Available in 10, 14 and 18 inch sizes. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Outer Layer	100% smooth knit polyester
Inner Layer	Urethane foam – polyester, yellow
Dimensions:	T6-10: 2.5 in X 10 in (6 cm x 25 cm) T6-14: 2.5 in X 14 in (6 cm x 35 cm) T6-18: 2.5 in X 18 in (6 cm x 45 cm)
For Use with: (as needed)	QDT TMop Adapters with SSU handles; QPT TMop Adapters with STA handles T6-18: QPWT Window Washing Tool with STA handles
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Packaging:	Individually bagged
Traceability:	From TAA: Lot number and part number printed on each bag and case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website