

ISOLATOR TOOL COVER, MEGATEX

Product Specification

Product Name:	Isolator Tool Cover, MegaTex
Manufacturer Part Number(s):	ITC620-9, ITC620-911, ITC620-9NS, ITC620-911NS
Description:	A slip-on non-woven polyamide mop cover to be used with QDIT or QPIT Isolator Tools. Manufactured at Micronova's Torrance, CA facility. Designed for use in biosafety cabinets, hoods, and other small clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Cleaning Surface	Non-woven polyamide
Lining	100% polyester
Attachment:	Elastic Band
Interior	NS suffix: No elastic ITC620-911 and ITC620-911NS: urethane foam ITC620-9 and ITC620-9NS: no foam
Dimensions:	2 in x 9 in (5 cm x 22.5 cm)
For Use with: (as needed)	QDIT Stainless Steel Isolator Tool with SSU series handles or QPIT Plastic Isolator Tool with STA series threaded handles.
Autoclavability:	Autoclavable
Chemical Compatibility:	Good Compatibility with standard disinfectants and cleaning solutions
Packaging:	Individually bagged
Traceability:	Part number and lot number on each bag and case
Irradiation Standards:	Not Irradiated
Certificates:	Certificate of Conformance available upon request and on website