



MEGATIP POLYESTER IRRADIATED

Product Specification

Product Name:	MegaTip Polyester
Manufacturer Part Numbers:	MT-6IR
Description:	A polyester/polyethylene foam slipcover with stainless steel fasteners for use with the QDMT-1 series mop adapter. Manufactured at Micronova's Torrance, CA facility Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001:2008 Facility
Physical Data:	
Outer layer:	100% Polyester
Inner layer:	Polyethylene Foam
Measurement:	
MT-6IR:	1-1/2" Diameter x 4" Long
Packaging:	Each mop head is packaged in a Class 100 clean room. Double bagged. Irradiation indicator placed on outer bag. 12 mops are then packed into a poly lined case. Irradiation indicator placed on case.
Traceability:	Lot numbers printed on each outer bag and each case
Irradiation Specification:	Minimum irradiation dose 40 kGy
Certificates:	Certificate of irradiation shipped with each order