

## CLEANROOM PEN

### Product Specification

Product Name:	Cleanroom Pen
Manufacturer Part Number(s):	PEN-01, PEN-02, PEN-03
Description:	A low-sodium ink, capped ballpoint pen designed with a special poly barrel to limit particle generation. Manufactured at Micronova's Torrance, CA facility. Designed for clean and controlled environments.
Manufacturing Standards:	ISO 9001 Certified Facility
Physical Data:	
Materials of Construction:	
Barrel:	White high-density poly barrel
Ink:	Low-sodium ink PEN-01: black ink PEN-02: blue ink PEN-03: red ink
Dimensions:	
Length x Width:	5¼ in x ¾ in
Nib:	0.7 mm
For Use with: (as needed)	Cleanroom documentation
Autoclavability:	Not autoclavable.
Chemical Compatibility:	Holds up to wipe down with IPA, bleach and phenols. Integrity of printed barrel will be compromised with thinners and strong solvents.
Packaging:	Packaged 10 pens per bag, and then double bagged.
Traceability:	Lot numbers printed on each outer bag and each case.
Certificates:	Certificate of Conformance available upon request.