



## CLEANROOM PEN - IRRADIATED

### Product Specification

Product Name:	Cleanroom Pen - Irradiated
Manufacturer Part Numbers:	PEN-01TIR, PEN-02TIR, PEN-03TIR
Physical Data:	White high-density poly barrel Low-sodium ink Available in black, red and blue
Chemical Resistance:	Holds up to wipedown with IPA, bleach and phenols. Integrity of printed barrel will be compromised with thinners and strong solvents.
Manufacturing Standards:	ISO 9001:2008 facility
Packaging:	Each pen individually bagged. Pens are then triple bagged in groups of 10. 100 pens per poly lined case.
Traceability:	Lot numbers printed on each outer bag and each case
Irradiation Specification:	Minimum irradiation dose 40 kGy
Certificates:	Certificate of irradiation shipped with each order

***Pens will hold up to wipedown with IPA, bleach and phenols.  
Integrity of printed barrel will be compromised with thinners and strong solvents.***